

# Amplisure<sup>®</sup> HCV Quantitative Kit (Real Time PCR Kit)

**IVD**

**Product Insert**



**RAS Lifesciences Pvt. Ltd**

Plot No. 13, 4-7-18/13/2., Raghavendra Nagar,

Nacharam, Hyderabad 500 076, India

Tel: +91-40-65261562,

[www.raslifesciences.com](http://www.raslifesciences.com)

## Index

	<b>Page No.</b>
Introduction	4
Product Description	4
Recommended Work Areas	4
General Precautions	5
Precautions while extracting Nucleic acid	5
Precautions while setting up a PCR reaction	5
Precautions for post PCR or equipment area/room	5
Precautions after completion of Real time PCR assay	5
Usage Limitations	6
Safety Precautions	6
Storage Conditions and Product Stability	6
Symbols	7
Kit Components	8
Materials required but not supplied	8
Quality Systems	9
Sample Type/Collection/Storage/Transport	9
Sample Type	9
Sample Collection, Storage and Transport	9
<b>Assay Procedure</b>	<b>10</b>
RNA Extraction	10
Use of Internal Control (IC)	10
Usage of Internal Control at the RNA extraction step	11
Usage of Internal Control at the Real time PCR step	11
qPCR Protocol	11
Preparation of Reaction Mastermix	11
PCR Programming	13
Plate Setup	13
Cycling conditions	14
Data Analysis	14
Setting the threshold for the qPCR Data analysis	15
Results	15
Interpretation	15
Troubleshoot	16

<b>Assay Characteristics</b>	<b>17</b>
Analytical Sensitivity	17
Analytical Sensitivity of qPCR assay	17
Analytical Sensitivity in consideration of Extraction Method	17
Linear Range	18
Conversion Factor (IU/ $\mu$ L to copies/ $\mu$ L)	18
Specificity	19
Detection of Different HCV genotypes (1-6) using International Standard	19
Cross Reactivity Data	20
Precision	21
Intra-assay variability	21
Inter-assay variability	22
<b>Abbreviations</b>	<b>23</b>
<b>References</b>	<b>24</b>
<b>Ordering Information</b>	<b>24</b>

## Introduction

Hepatitis C virus (HCV) is a hepatotropic virus of family Flaviviridae and genus Hepacivirus with a single strand positive RNA genome. It is a leading cause of chronic liver disease and has already infected at least 170 million people worldwide. In India, there are about 12-13 million HCV carriers and modeling data predict that the burden of disease could soon increase substantially. Till date, six different HCV genotypes and more than 70 subtypes have been identified based on the nucleic acid sequences. Although most chronic HCV patients have mild chronic hepatitis, it is a progressive disease and can lead to cirrhosis or hepatocellular carcinoma. A large number of genotypes have been identified among hepatitis C virus isolates from all over the world. It has been suggested that different genotypes have different clinical outcomes with regard to disease severity and response to interferon therapy.

## Product Description

Amplisure<sup>®</sup> HCV Quantitative Kit is an *in-vitro* diagnostic kit for quantitation of Hepatitis C Virus (1-6 genotypes) in human plasma. The kit contains the necessary reagents for performing HCV quantitation by Real-time PCR.

Pathogen detection by Real Time polymerase chain reaction (PCR) is based on the amplification of specific region (5'UTR) of the pathogen genome. The assay principle is based on Taqman probes which allow higher specificity and sensitivity.

In addition, the Amplisure<sup>®</sup> HCV Quantitative Kit contains a second amplification system to identify possible PCR inhibition by using an internal control (IC), without affecting the analytical sensitivity of the assay. External quantitation standards calibrated with WHO control are supplied, which allow the determination of the amount of Hepatitis C Viral RNA.

## Recommended Work areas

Molecular Diagnostics work area includes:

- a) Sample preparation area/room – for extraction of nucleic acids from clinical samples
- b) Pre-PCR area/room - for setting up PCR reaction
- c) PCR area/room – for performing PCR using the thermocyclers

As part of Good Laboratory Practices (GLP), it is recommended to have dedicated areas to avoid cross contamination.

## General Precautions

### ***Precautions while extracting Nucleic acid***

Always wear proper attire (powder free gloves, facemask and Head cap) before starting the nucleic acid extraction procedure. During preparation of samples, compliance with good laboratory practices are essential to minimize the risk of cross-contamination between samples, and the inadvertent introduction of ribonucleases (RNases) into samples during and after the extraction procedure. Proper aseptic technique should always be used when working with RNA.

The Sample Preparation Area is dedicated to processing samples. All reagents used in the Sample Preparation Area should remain in this dedicated area at all times. Laboratory coats, pipettes, pipette tips and vortex mixer used in the Sample Preparation Area must remain in this area and not be moved to the Pre-PCR/PCR area. Discard the gloves before leaving this area. Do not bring amplified product into the Sample Preparation Area. Usage of filter tips is recommended while sample preparation and should be performed in a Biosafety cabinet.

### ***Precautions while setting up a PCR reaction***

PCR assay is sensitive and any accidental introduction of product from previous amplification reactions leads to incorrect results. Hence measures to reduce the risk of contamination in the laboratory should include physically separating the activities involved in performing PCR and complying with good laboratory practices.

It is recommended to have proper cleaning procedures to minimize the risk of cross contamination and carry over contamination (e.g. DNA OUT™, RNase OUT™, 0.1% Sodium Hypochlorite, Fumigation etc.).

It is recommended that areas should be defined in Pre-PCR room for preparation of mastermix and addition of templates. Laboratory coats and equipment used in the Pre-PCR Area must remain in this area and should not be moved to the Sample Preparation Area.

### ***Precautions for post PCR or equipment area/room***

The Real time PCR instrument/s should be kept in a separate segregated area away from Sample preparation area and Pre-PCR area.

### ***Precautions after completion of Real time PCR assay***

The reaction tubes or strips should be properly discarded without opening the caps, after the completion of run to avoid carry over contamination.

## Usage Limitations

1. The kit and all its components are for *in-vitro* diagnostics only.
2. The product is to be used by personnel specially trained in the *in-vitro* diagnostics procedures only.
3. Follow the product insert strictly for optimal PCR results.
4. Do not use the kit beyond the expiry date mentioned on the kit box.
5. Follow the guidelines provided in product insert for sample collection, storage and transport.
6. For ideal performance, store the kit under recommended conditions only.











## Safety Precautions

1. All patient specimens should be considered as potentially infectious and handled in a BSL2 biosafety hood with BSL3 practices.
2. Wear personal protective equipment, including gloves, head cap, face mask and lab coats when handling kit reagents/sample extraction. Wash hands thoroughly using detergents before and after performing the test.
3. Do not smoke, drink or eat in areas where kit reagents and/or human specimens are being used.
4. Dispose of unused kit reagents and human specimens as per regulatory guidelines.

## Storage Conditions and Product Stability

1. All the kit reagents should be stored at -20°C. Replace all the kit components immediately at -20°C after usage.
2. Repeated thawing and freezing (more than 6 x) of all kit reagents should be avoided, as it reduces assay sensitivity. If needed, make aliquots of the kit reagents according to the volume used in the protocol prior to freezing.
3. Allow reagents to be thawed completely on Ice/4°C prior to use.
4. Kit reagents are stable through the end of the expiration date indicated on the box when stored at -20°C.

## Symbols

Description of Symbol	Denotation
	<i>in-vitro</i> Diagnostic medical device
	Consult Instruction manual (Product Insert) for use
	Lot Number of the kit or Kit contents
	Catalogue number of Kit
	Contains sufficient for <N> reactions (Pack Size)
	Manufacturer
	Temperature limitation (Storage Condition)
	Use by MMM-YYYY (Expiry Date)
	Biological risk (handle carefully)
	Important Note

## Kit Components

Color Coding (Caps)	Contents	Description	50 rxns (QT-HCV-50)	100 rxns (QT-HCV-100) 2 x 50 rxns
Yellow	RAS qRNA PCR Mix	DNA Amplification Reagent	1 x 750 µL	2 x 750 µL
Red	RAS RT Mix	cDNA synthesis Reagent	1 x 350 µL	2 x 350 µL
Brown	RAS HCV PPM	HCV Primer-Probe Mix	1 x 110 µL	2 x 110 µL
Lavender	RAS HCQS1 (2 x 10 <sup>4</sup> IU/ µL)	HCV Quantitation Standards	1 x 60 µL	2 x 60 µL
Lavender	RAS HCQS2 (2 x 10 <sup>3</sup> IU/ µL)		1 x 60 µL	2 x 60 µL
Lavender	RAS HCQS3 (2 x 10 <sup>2</sup> IU/ µL)		1 x 60 µL	2 x 60 µL
Lavender	RAS HCQS4 (2 x 10 <sup>1</sup> IU/ µL)		1 x 60 µL	2 x 60 µL
Natural	RAS IC PCR Mix	Internal Control	1 x 50 µL	2 x 50 µL
Natural	RAS R-IC EX Mix		1 x 500 µL	2 x 500 µL
White	MBGW	Molecular Biology Grade water	1 x 1.0 mL	2 x 1.0 mL

## Materials required but not supplied

The materials which are required but not supplied are listed below:

1. Viral RNA Extraction kit
2. Biosafety Cabinet
3. PCR Hood
4. Calibrated variable micropipettes
5. Sterile pipette filter tips (aerosol free)
6. Vortex mixer
7. Dry Bath
8. Benchtop centrifuge with rotor for 1.5 mL reaction tubes
9. Real Time PCR machine



10. Strip Tubes and Caps (0.2 mL) or PCR Tubes (0.2 mL) or 96 well plate
11. Cooling block (96 x 0.2 mL tubes)
12. 1.5 mL centrifuge tubes
13. 1.5 mL centrifuge tube stand
14. Cooling block (2 mL tubes)
15. Sterile powder free gloves
16. Facemask
17. Head cap
18. Lab coats

## Quality Systems

In accordance with ISO-certified Quality Management System (9001:2008 and 13485: 2003) of RAS Lifesciences, each lot of Amplisure<sup>®</sup> HCV Quantitative PCR Kit is tested against predetermined specifications to ensure consistent product quality.

## Sample Type/Collection/Storage/Transport

### Sample Type

Plasma (K<sub>2</sub>EDTA-Blood)

*Heparinized Blood must not be used as they inhibit the PCR reaction*

### Sample Collection, Storage and Transport

3-5 mL of blood has to be drawn into a K<sub>2</sub>EDTA vacutainer. Cap and swirl the tubes for uniform mixing of blood and K<sub>2</sub>EDTA. When K<sub>2</sub>EDTA is used, whole blood can be collected in tubes with or without a gel separator. The K<sub>2</sub>EDTA blood samples are centrifuged (20 minutes at 1000–1500 x *g*) to separate plasma from cellular material and in the case of a non gel separator tube, the plasma should be removed to a secondary sterile tube within 4 hours of phlebotomy.

Collection and storage of unstabilized whole blood is not recommended for PCR analysis, because RNA degradation occurs in blood stored *ex vivo*. The sensitivity of the assay can be reduced if whole blood samples are frozen or stored for a longer period of time.

Plasma separated in a gel separator tube may be transported to the laboratory *in situ*. Plasma should be shipped at 2 to 8 °C and stored at -20°C as it is stable for up to five days at 2 to 8 °C and longer if frozen at -20°C or -70°C or lower. Do not store plasma samples in a “frost -free” freezer as the temperature is cycled several times per day on this type of freezer, causing degradation of nucleic acid targets.

Sample material should be transported in a leak proof, unbreakable transport container to avoid leakage of sample. The samples should be transported following the local and national instructions for the transport of pathogen material.

## Assay Procedure

### RNA Extraction

Amplisure<sup>®</sup> HCV Quantitative Kit has been validated using the following Viral RNA extraction kits:

- 1- Roche High Pure Viral RNA kit (Cat. No. 11858882001)
- 2- QIAamp Viral RNA Mini Kit (Cat. No. 52904)

Follow the manufacturer's instructions mentioned in the manual for Viral RNA extraction. Different pack sizes of the above mentioned kits can be used. However the customer can also validate their own extraction process using other Viral RNA extraction Kits.

Recommended sample volume for extraction and elution are as follows:

Sl. No.	Name of the RNA Isolation Kit	Recommended Sample volume ( <i>to be taken for RNA Extraction</i> )	Recommended Final Elution volume
1.	Roche High Pure Viral RNA kit (Cat. No.11858882001)	200 µL	50 µL
2.	QIAamp Viral RNA Mini Kit (Cat. No. 52904)	140 µL	50 µL

The analytical sensitivity of the assay in consideration of the purification was determined using the above defined volumes.

### Use of Internal Control (IC)


Internal controls are supplied (RAS R-IC Ex Mix and RAS IC PCR Mix) along with Amplisure<sup>®</sup> HCV Quantitative Kit. This allows the user to control the Viral RNA isolation procedure as well as to check for possible PCR inhibition.



Internal control should only be used once, either at the Extraction step or at the PCR step

### ***Usage of Internal Control at the RNA extraction step***

If internal control (IC) is required to be added at the time of RNA extraction, add 10 µL of RAS *R-IC Ex mix* per isolation to the lysis buffer along with other components of kit used for lysis (as per kit instructions) and vortex for 5 seconds prior to usage.



**Fig. 1- Viral RNA Extraction Overview**

### ***Usage of Internal Control at the Real time PCR step***

The internal control can optionally be used exclusively to check for possible PCR inhibition. For this application, add the internal control directly to the PCR master mix as described on Pg. No 15. Addition of Internal control should be done only once, either at the time of sample extraction or during PCR setup.

## **qPCR Protocol**

### **Preparation of Reaction Mastermix**

Quantitation procedure with Amplisure<sup>®</sup> HCV Quantitative Kit involves *1 step RT qPCR*. It is recommended that a minimum of three standards and a negative control (MBGW should be used as negative control) are required to be included in a single run for acquiring proper results.

Set up a real time single step RTPCR reaction as follows:

**1- qPCR reaction mix composition without Internal Control (When R-IC Ex Mix is added during RNA extraction)**

<b>Components</b>	<b>Volume per reaction ( µL) (for final vol. of 30 µL)</b>
RAS qRNA PCR Mix	15.0
RAS RT Mix	7.0
RAS HCV PPM	2.0
RNA/HCQS/ MBGW	5.0
MBGW	1.0

**2- qPCR reaction mix composition with Internal Control (When R-IC Ex Mix is not added during RNA extraction)**

<b>Components</b>	<b>Volume per reaction ( µL) (for final vol. of 30 µL)</b>
RAS qRNA PCR Mix	15.0
RAS RT Mix	7.0
RAS HCV PPM	2.0
RNA/HCQS/ MBGW	5.0
RAS IC PCR Mix	1.0



- 1. Addition of IC PCR mix (if required) should be done at the time of preparation of master mix.**
- 2. The results may be inconsistent, if the IC PCR mix is added individually.**

Place the PCR plate/tubes/strips in real time thermocycler.

## PCR Programming

The Amplisure HCV Quantitation kit is validated on the following instruments:

- Rotor-Gene™ 6000
- Rotor-Gene™ Q 5plex
- ABI 7500 DX Real-Time PCR System
- ABI 7300 Real-Time PCR System
- Eppendorf Realplex 4
- Bio-Rad™ CFX 96

### Plate Setup

1. Program the plate setup by labeling the slots as per tube/strip/plate labels. The sequence of labeling of slots should be the same way as the tube/strip/plate is kept in the machine.
2. Select the type of sample (Unknown/Standard/NTC) for each slot.
3. When Standard is selected as sample type, mention the quantitation value in IU/μL (2E4, 2E3 etc)
4. Select the channel for acquisition (FAM/Yakima Yellow)

Sl. No.	Name of channel	Source wavelength (nm)	Detection wavelength (nm)
1.	FAM (Pathogen target)	470	510
2.	Yakima Yellow (Internal Control)	530	555

5. For background calibration in different instruments, follow the procedure described below:
 

Rotor-Gene™ 6000	- Perform 'Gain optimization'
Rotor-Gene™ Q 5plex	- Perform 'Gain optimization'
ABI 7500 DX Real-Time PCR System	- Select Passive Reference dye 'ROX'
ABI 7300 Real-Time PCR System	- Select Passive Reference dye 'ROX'
Eppendorf Realplex 4	- Select 'ROX' for background calibration
Bio-Rad™ CFX 96	- Select 'ROX' for background calibration



**Preparation of reaction mastermix and cycling conditions are same for all the instruments listed in the product insert. For instrument specific protocols, please contact our technical support team at [amplisure@raslifesciences.com](mailto:amplisure@raslifesciences.com)**

### ***Cycling conditions***

1. Configure the following program in the machine.

<b>Steps</b>	<b>No. of cycles</b>	<b>Temperature (°C)</b>	<b>Time</b>
1 (Reverse Transcription)	1	42	15 min.
2 (Initial denaturation)	1	95	10 min.
3 (PCR cycling)	40	95	15 sec.
		60	30 sec
		72*	30 sec
<b>* Plate Read/Data Acquisition in FAM and Yakima Yellow channel</b>			

2. Set the reaction volume as 30 µL.
3. Plate read/Data Acquisition for FAM and Yakima Yellow channel should be incorporated in the third stage of step 3 (72°C/30 sec).
4. The ideal run time for the assay is 120 minutes. Note: *In case of Eppendorf Realplex 4, select RAMP rate as 35%.*



**Preparation of reaction mastermix and cycling conditions are same for all the instruments listed in the product insert. For instrument specific protocols, please contact our technical support team at [amplisure@raslifesciences.com](mailto:amplisure@raslifesciences.com)**

### **Data Analysis**

Analyze the data after completion of the run. Check the  $R_n/Cycle$  amplification plot and  $R_n/Cycle$  amplification plot to observe the amplification signal generated by different samples in the run. Compare both the plots for data analysis. Also look for noisy signals, if observed as it might not give you a proper result.

## Setting the threshold for the qPCR Data analysis

The threshold should be set either automatically (by the machine itself)/ or manually just above the background signal of the negative controls and negative samples by referring to  $R_n$ /Cycle amplification plot. The mean threshold value calculated from these experiments will most likely work for the majority of future runs, but the user should nevertheless review the generated threshold value at regular intervals.

## Result

The values for unknown samples would appear in the result column in  $IU/\mu L$  in FAM Channel. Samples showing no amplification in FAM channel should show amplification in VIC channel, and then only results should be considered. The negative control should not show any value in the result column.

## Interpretation

Interpret the values for unknown samples, only if the Slope of Standards is between -3.1 to -3.6 and PCR efficiency is between 90%-110% (0.9 - 1.1) and there should be no amplification in negative control.

Observation	Interpretation	Conclusion
Amplification signal detected in HCV channel (FAM) and in Internal control channel (Yakima Yellow)	HCV RNA within quantitation range	Proceed for further Analysis
Amplification signal detected in HCV channel (FAM) but no signal in Internal Control channel (Yakima Yellow)	HCV RNA within quantitation range	
Amplification signal not detected in HCV channel (FAM) but detected in Internal control channel (Yakima Yellow)	HCV RNA below quantitation limit	
No Amplification signal detected in HCV channel (FAM) as well as Internal control channel (Yakima Yellow) in unknown samples	Possible inhibition of PCR	Dilute the RNA sample (1:100) and repeat the assay*

To convert the results from IU/μL to IU/mL use the following formula:

$$\text{Result (IU/mL)} = \frac{\text{Result (IU/}\mu\text{L)} \times \text{Elution Volume (}\mu\text{L)}}{\text{Sample Volume (mL)}}$$

**i** \* For calculating the result of diluted sample (1:100); multiply the observed IU/mL value by dilution factor, 100

$$\text{Result of 1:100 diluted sample (IU/mL)} = \text{Dilution Factor x } \frac{\text{Result (IU/}\mu\text{L)} \times \text{Elution Volume (}\mu\text{L)}}{\text{Sample Volume (mL)}}$$

### Troubleshoot

Observation	Possible cause	Solution
No amplification signal for Standards in FAM channel	<ol style="list-style-type: none"> <li>One of the components may not have been added.</li> <li>Incorrect channel selection</li> <li>Incorrect programming of the real time machine.</li> <li>Instrument is not working properly</li> </ol>	<ol style="list-style-type: none"> <li>Repeat the assay by following the correct protocol and addition of reagents</li> <li>2&amp;3. Please recheck the PCR program</li> <li>4. Contact manufacturer of thermocycler for technical support.</li> </ol>
Weak amplification signal for standards (Signal below threshold) in FAM Channel	<ol style="list-style-type: none"> <li>Improper PCR programming.</li> <li>Inaccurate dispensing of reagents</li> <li>Possible deterioration of kit components due to improper storage</li> </ol>	<ol style="list-style-type: none"> <li>Repeat the assay by following the correct protocol</li> <li>Minimize Pipetting errors/Check for calibration status of pipettes</li> </ol>
Identical/Similar Ct values observed in FAM channel	<ol style="list-style-type: none"> <li>Possible contamination of Kit reagents / Standards/Work area.</li> </ol>	<ol style="list-style-type: none"> <li>Use fresh aliquots of Standards/Kit Reagents (if available)</li> <li>Clean the PCR rack/Pipettes thoroughly as per GLP</li> <li>Clean and Fumigate the work area overnight prior to use</li> </ol>

**i** For any other technical query; please contact [amplisure@raslifesciences.com](mailto:amplisure@raslifesciences.com)



## Assay Characteristics

### Analytical Sensitivity

The analytical sensitivity of PCR assay as well as the analytical sensitivity in consideration of the purification (RNA Extraction) was assessed for the Amplisure<sup>®</sup> HCV Quantitative Kit.

The analytical sensitivity in consideration of the purification is determined using WHO standards/HCV positive clinical samples in combination with a particular extraction method.

### Analytical Sensitivity of qPCR assay

The analytical sensitivity is determined independent from the selected extraction method, using an international standard of known concentration ie. 4<sup>th</sup> WHO International Standard for Hepatitis C Virus for Nucleic acid amplification techniques (NIBSC code: 06/102).

To determine the analytical sensitivity of the Amplisure<sup>®</sup> HCV Quantitative Kit, a standard dilution series was set up from 1.04 IU/ $\mu$ L to 41.6 IU/ $\mu$ L and analyzed with Amplisure HCV Real Time PCR Quantitation Kit.

Testing was carried out on 10 different days on 3 replicates. The results were analyzed by statistical analysis.

***The analytical sensitivity of the Amplisure HCV Quantitative Kit is 2.08 IU/ $\mu$ L.***

### Analytical Sensitivity in consideration of Extraction Method

The analytical sensitivity in consideration of the purification (Roche High Pure Viral RNA kit; Cat. No.11858882001 and QIAamp Viral RNA Mini Kit; Cat. No. 52904) of the Amplisure HCV Quantitative Kit was also determined using a dilution series of the 4<sup>th</sup> WHO International Standard for Hepatitis C Virus for Nucleic acid amplification techniques (NIBSC code: 06/102) from 2.6 to 2600 IU/mL spiked in healthy plasma specimens (Negative for HCV by Real time PCR).

These were subjected to RNA extraction using the Roche High Pure Viral RNA kit and QIAamp Viral RNA Mini Kit (Roche, extraction volume: 0.2 mL, elution volume: 50  $\mu$ L and Qiagen, extraction volume: 0.14 mL, elution volume: 50  $\mu$ L)

Each of the dilutions was analyzed with the Amplisure<sup>®</sup> HCV Quantitative Kit on 10 different days on 3 replicates. The results were determined by statistical analysis.

***The analytical sensitivity in consideration of the purification of the Amplisure<sup>®</sup> HCV Quantitative Kit, is 260 IU/mL (with Roche High Pure Viral RNA kit) and is 300 IU/mL (with QIAamp Viral RNA Mini Kit)***

## **Linear Range**

The linear range (analytical measurement) of the Amplisure<sup>®</sup> HCV Quantitative Kit was determined by analyzing a dilution series of HCV quantitation standards from  $2 \times 10^7$  IU/ $\mu$ L to 2 IU/ $\mu$ L. The dilution series has been calibrated against the 4<sup>th</sup> WHO International Standard for Hepatitis C Virus for Nucleic acid amplification techniques (NIBSC code: 06/102).

Each dilution was tested in replicates (n= 3) using the Amplisure<sup>®</sup> HCV Quantitative Kit.

The linear range of the Amplisure HCV Real Time PCR Quantitation Kit has been determined to cover concentrations from **20 IU/ $\mu$ L to  $2 \times 10^7$  IU/ $\mu$ L**.


## **Conversion Factor (IU/ $\mu$ L to Copies/ $\mu$ L)**

The 4<sup>th</sup> WHO International Standard for Hepatitis C Virus for Nucleic acid amplification techniques (NIBSC code: 06/102) has been assigned an International Unit (IU) by the WHO as the number of genome equivalents.

The conversion factor value for Amplisure<sup>®</sup> HCV Quantitative kit has been assigned after calibration of the HCQS (HCV quantitation standards) with the 4<sup>th</sup> WHO International Standard for Hepatitis C Virus for Nucleic acid amplification techniques (NIBSC code: 06/102). The Conversion Factor (IU/ $\mu$ L to Copies/ $\mu$ L) is entirely dependent on the current PCR assay kit. A Correlation and Regression analysis was performed to compare the two standard curves of HCV RNA values obtained by Amplisure HCV Quantitative kit (copies/ $\mu$ L) and HCV NIBSC standards (IU/ $\mu$ L).

The correlation coefficient between the expected and estimated values was very good  $r = 0.9912$ . Please see Fig. 2 below

$r = 0.9912$



**Fig2. Correlation and Regression plot of HCV NIBSC standards (IU/μL) and HCV RNA (copies/μL)**

The conversion factor calculated for Amplisure HCV PCR kit from IU/μL to Copies/ μL is 2.18, i.e. **1 IU/μL = 2.18 copies/ μL**

To convert into copies/mL use the following formula:

$$\text{Result (copies/mL)} = \frac{\text{Calculated Value (copies/}\mu\text{L)} \times \text{Elution Volume (}\mu\text{L)}}{\text{Sample Volume (mL)}}$$


### Specificity

The specificity of the Amplisure<sup>®</sup> HCV Quantitative Kit is ensured by the selection of the primers and probes, as well as the selection of stringent reaction conditions. The primers and probes were checked for possible homologies to all published sequences (Genbank) by BLAST analysis.

## Detection of Different HCV genotypes (1-6) using International Standard

The detectability of all relevant genotypes (1-6) was ensured by a Database alignment and by a qPCR run with the following genotypes, using International standard (2<sup>nd</sup> HCV RNA genotype panel for Nucleic Acid Amplification Techniques –NIBSC code: 08/264) and using ≥ 300 positive HCV samples covering all genotypes.

Sl. No.	International HCV Genotype Controls	Source of International controls	Result obtained with Amplisure <sup>®</sup> HCV Quantitative Kit
1.	Geno 1	NIBSC	Detected
2.	Geno 2	NIBSC	Detected
3.	Geno 3	NIBSC	Detected
4.	Geno 4	NIBSC	Detected
5.	Geno 5	NIBSC	Detected
6.	Geno 6	NIBSC	Detected



**Fig3: Detection of HCV RNA Genotype panel (NIBSC) using Amplisure HCV Quantitative Kit**

## Cross Reactivity Data

A potential cross-reactivity of the Amplisure HCV Quantitative Kit was tested using the control group listed below. None of the tested pathogens has been reactive. No cross-reactivities appeared with mixed infections.

The specificity was also validated with 50 different healthy plasma specimens/other various sample types. These did not generate any amplification signals with the HCV specific primers and probes, which are included in the Amplisure HCV Quantitative Kit.

<b>Pathogen Tested</b>	<b>Cross Reactivity with the HCV Primers/Probes</b>
Adenovirus	-ve
Herpes Simplex Virus-1	-ve
Herpes Simplex Virus-2	-ve
Epstein Barr virus	-ve
Human Immunodeficiency Virus	-ve
Cytomegalo virus	-ve
Hepatitis B virus	-ve
Enterovirus	-ve
BK virus	-ve
MTB complex	-ve
Plasmodium sps.	-ve

## **Precision**

The precision data of the Amplisure<sup>®</sup> HCV Quantitative Kit have been generated for HCV positive clinical specimens and HCV quantitation standards.

### **Intra-assay variability**

The precision data consists of the *Intra-assay variability* which is variability of multiple results of samples of the same concentration within one experiment.

The data obtained were used to determine the standard deviation and the coefficient of variation for the pathogen specific PCR. Precision data of the Amplisure<sup>®</sup> HCV Quantitative Kit have been collected using all the quantitation standards (HCQS1- HCQS4). Also included were three different samples with different viral loads. Testing was performed with 3 replicates for 10 days.

## Inter-assay variability

The precision data consists of the *Inter-assay variability* which is variability of multiple results of samples of the same concentration within different experiments.

The data obtained were used to determine the standard deviation and the coefficient of variation for the pathogen specific PCR. Precision data of the Amplisure® HCV Quantitative Kit have been collected using all the quantitation standards (HCQS1- HCQS4). Also included were three different samples with different viral loads. Testing was performed with 3 replicates for 10 days.

The precision data was calculated on basis of the CT values obtained.

Sample Type	Mean [Log Value (IU/μL)]	Variability Testing	Standard Deviation	Coefficient of Variation (%)
HCQS1	Log 4.3	Intra-assay variability	0.012	0.001
		Inter-assay variability	0.090	0.004
HCQS2	Log 3.3	Intra-assay variability	0.009	0.001
		Inter-assay variability	0.080	0.003
HCQS3	Log 2.3	Intra-assay variability	0.009	0.001
		Inter-assay variability	0.250	0.009
HCQS4	Log 1.3	Intra-assay variability	0.020	0.001
		Inter-assay variability	0.050	0.001
Control 1	Log 5.2	Intra-assay variability	0.190	0.008
		Inter-assay variability	0.290	0.012
Control 2	Log 3.7	Intra-assay variability	0.100	0.003
		Inter-assay variability	0.220	0.007
Control 3	Log 2.2	Intra-assay variability	0.140	0.004
		Inter-assay variability	0.340	0.010

## Abbreviations

<b><i>Abbreviation</i></b>	<b><i>Expansion</i></b>
HCV	Hepatitis C Virus
RNA	Ribonucleic Acid
IC	Internal Control
RNases	Ribonucleases
PCR	Polymerase Chain Reaction
BSL2	Bio Safety Level 2
BSL3	Bio Safety Level 3
mL	Milli Liters
µL	Micro Liters
K <sub>2</sub> EDTA	Potassium Ethylene Diamine Tetra Acetate
g	Relative Centrifugal Force
qPCR Protocol	Quantitative PCR protocol
MBGW	Molecular Biology Grade Water
RT PCR	Real Time PCR
NTC	No Template Control
FAM	Carboxyfluorescein
ROX	Carboxy-X-rhodamine
NIBSC	National Institute for Biological Standards and Control
IU	International Units
WHO	World Health Organization
BLAST	Basic Local Alignment Search Tool
MTB complex	Mycobacterium tuberculosis complex
Sps	Species
Rxn	Reaction

## References

1. Handbook "Sampling, Transportation, and Storage of Clinical Material for PCR Diagnostics", developed by Federal State Institution of Science Central Research Institute of Epidemiology of Federal Service for Surveillance on Consumers' Rights Protection and Human Well-Being, Moscow, 2008.
2. Machnik G, Pelc E, Zapala M, Gasecka-Czapla M, Kaczmarczyk G, Kozłowska D, Okopien B. Designing and optimization of real-time RT-PCR technique for the detection of hepatitis C virus (HCV) genome in blood serum as internal laboratory quality control. *Przegl Epidemiol.* 2011; 65(2):325-32.
3. Narahari S, Juwle A, Basak S, Saranath D. Prevalence and geographic distribution of Hepatitis C Virus genotypes in Indian patient cohort *Infect Genet Evol* 2009;(4):643-5

## Ordering Information



QT-HCV-50	: 50 rxns
QT-HCV-100	: 100 rxns