

RAS Real *Taq* DNA Polymerase

(Modified DNA Polymerase)

5,000 U/ml

1600 units

Lot no xxxxxxxxx

Exp bcbcbcc

Store at -20C

Description

RAS Real *Taq* DNA Polymerase is a modified form of thermostable DNA polymerase from the thermophilic bacterium *Thermus aquaticus*. The enzyme catalyzes 5'→3' synthesis of DNA, has no detectable 3'→5' exonuclease (proofreading) activity and possesses 5'→3' exonuclease activity. It possesses high specificity and subsequently increases the yield of the specific PCR product. This is due to the reduced mis-primed products and primer dimers at ambient temperature as it is provided in an inactive state, thus it requires incubation at 95 °C for 15 minutes for activation.

Source

Produced from an *E. coli* strain that possesses the *Taq* DNA Polymerase gene from *Thermus aquaticus* that subsequently modified chemically.

Applications

- Routine PCR amplification of DNA fragments
- Quantitative PCR
- Nested PCR
- RFLP

Supplied in: 20mM HEPES, 100mM KCl , 0.1 mM EDTA,0.5mM PMSF , 1 mM DTT, 50% Glycerol

Concentration

5U/μL

Reagents Supplied with Enzyme

- 50 mM MgCl₂
- 10x Hot Start Buffer
 - 200 mM Tris HCl (pH 8.8),
 - 500 mM KCl,
 - 160 mM (NH₄)₂SO₄
 - Other additives

Storage conditions

- The recommended storage condition is -20°C.
- Avoid keeping at room temperature for long and repeated freeze thaws.

Quality Control

Quality Control assays

Mis-priming assay

No detectable non-specific bands were observed in agarose gel after PCR carried out to amplify E.coli gene from plasmid in Real *Taq* buffer using 2.5 units of Real *Taq* DNA polymerase and gene specific primers in presence of back ground human genomic DNA (250 ng/μl).

Amplification efficiency

Amplification efficiency (AE) of RAS Real *Taq* tested on real-time PCR using quantitative standards belonging to various viruses yielded a slope between -3.1 and -3.6 (AE= 0.9-1.1).

DNase activity

No detectable change in banding pattern was observed in agarose gel after incubating the 25 units of Real *Taq* DNA Polymerase with 500 ng of phage DNA in RAS *Taq* reaction buffer at 37°C for 4 hrs.

RNase activity

Incubation of Human RNA (250 ng) with 25 units of Real *Taq* DNA polymerase in RAS *Taq* reaction buffer at 37°C for 1 hr results in no detectable change in banding pattern as determined by agarose gel electrophoresis.

E. coli DNA contamination

Taq polymerase is tested free of *E. coli* DNA contamination

Assay Conditions

Recommended PCR Reaction Mix (Real-Time)

Component	Volume (μL)
10x Hot Start reaction Buffer	2.5
MgCl ₂ (50 mM)	1.0
Primer-Probe Mix	2.0
Template DNA	variable
Hot Start <i>Taq</i> DNA Polymerase	0.5
Molecular Grade Water	Make up to 25

Recommended Thermocycling conditions

Operation	Temp	Time	Cycles
Initial denaturation **	95°C	15 min	1
Denaturation	95°C	15-30 s	35-45 cycles
Annealing& extension	Tm-5°C	60 s	

** Initial denaturation for 15 min at 95°C is necessary for activation of enzyme.

Recommended PCR Reaction Mix (End point)

RAS Real *Taq* DNA Polymerase

(Modified DNA Polymerase)

5,000 U/ml

1600 units

Lot no xxxxxxxxx

Exp bcbcbbc

Store at -20C

Component	Volume (μ L)
10x Hot Start reaction Buffer	2.5
MgCl ₂ (50 mM)	1.0
Forward primer (10 pmol/ μ L)	0.4-0.7
Reverse primer (10 pmol/ μ L)	0.4-0.7
Template DNA	variable
Hot Start <i>Taq</i> DNA Polymerase	0.5
Molecular Grade Water	Make up to 25

Recommended Thermocycling conditions

Operation	Temp	Time	Cycles
Initial denaturation **	95°C	15 min	1
Denaturation	95°C	30 s	30-35 cycles
Annealing	Tm-°C	30-60 s	
Extension	72°C	1 min/kb	
Final Extension	72°C	5-10 min	--

** Initial denaturation for 15 min at 95°C is necessary for activation of enzyme.

References

1. Chien, A., Edgar, D.B. and Trela, J.M. (1976) *J. Bact.*, 127, 1550-1557.
2. Kaledin, A.S., Sliusarenko, A.G. and Gorodetskii, S.I. (1980) *Biokhimiya*, 45, 644-651.
3. Lawyer, F.C. et al. (1993) *PCR Methods and Appl.*, 2, 275-287.
4. Longley, M.J., Bennett, S.E. and Mosbaugh D.W. (1990) *Nucleic Acids Res.*, 18, 7317-7322.